

## Primary Adult Naïve CD4+ CD45RA+ Cells

Prepared by: David Randolph (drdrdr@uab.edu) at University of Alabama, Birmingham

Goal: To obtain large numbers of highly pure primary CD4+ CD45RO- CD25- cells from adult peripheral blood, and prepare them as naïve, Th0, and Th1 cells. This protocol combines two separate magnetic bead-based purification strategies, a positive selection step with Dynal CD4 beads followed by a negative selection step with EasySep Naïve Human T CD4 T cell beads, to give desired purity but with higher yields than can be obtained with FACS sorting and cells are “untouched” with no antibodies bound to their surface.

### Reagents:

- 1) Dynal CD4 Positive Isolation Kit (Invitrogen)
- 2) EasySep Negative Selection Human Naïve CD4 T Cell Enrichment Kit (Stemcell)
- 3) 1 buffy coat (Research Blood Components)
- 4) Anti-CD25-biotin
- 5) RPMI + 1% FCS
- 6) T Cell Media = RPMI 1640 + 10%FCS + Pen/Strep + 2mM L-glut + 10mM Hepes + 50 $\mu$ M 2- $\beta$ ME
- 7) 2X Freezing media = 8ml FCS + 2ml DMSO
- 8) Human recombinant IL-12 (R&D Systems)

### Step 1: Positive Selection of CD4 Cells with Dynal Beads

- 1) Buffy coat should contain 0.5-1 x 10<sup>9</sup> cells. Of these around 10% should be CD4 T cells, so there are about 0.5-1 x 10<sup>8</sup> CD4 T cells in the pack.
- 2) Dilute Buffy coat 1 part cells to 2 parts PBE
- 3) Wash Dynal beads once with PBE. Add beads at ratio of 5 beads to 1 cell = 5 x 10<sup>8</sup> beads = 1.2 ml of beads.
- 4) Incubate for 20 minutes at 4°C with gentle mixing.
- 5) Place tube on magnet for 2 minutes.
- 6) Discard supernatant.
- 7) Gently wash bead-bound cells 3 times with cold PBE.
- 8) Resuspend the bead-bound cells in 10 ml of RPMI/1% FCS.
- 9) Release cells by adding 400 ul of Detach-a-bead solution and incubate for 45 minutes at room temp with gentle rotation.
- 10) Place on magnet to remove beads.
- 11) Wash beads with RPMI/1% FCS x 2 to remove as many cells as possible.
- 12) Collect all cells and count. Cells should be highly pure CD4+ cells.

### Step 2: CD45RO and CD25 depletion

- 1) Wash cells x 2 in PBE to remove Detachabead.
- 2) Per million purified CD4 T cells, add 0.5 anti-CD25-biotin.
- 3) Incubate on ice x 30 minutes.
- 4) Wash once gently with PBE.

- 5) Resuspend in PBE at  $5 \times 10^7$  cells/ml
- 6) Add 50ul CD45RO-biotin (from kit) per ml of cell suspension.
- 7) Incubate at RT for 15'
- 8) Add 50ul Enrichment Cocktail per ml of cells
- 9) Incubate at RT for 10min.
- 10) Add 100ul StemCell nanoparticles and incubate for 10min
- 11) Place tube in magnet for 10 min.
- 12) Pour of unbound CD4 cells.
- 13) Wash beads and repeat.
- 14) Count final cells, spin down, and resuspend in T Cell media.

### Step 3: Stimulation for DNase I experiments

- 1) For Naïve Cells, take 10 million in T Cell Media. Add 2X Freezing media and freeze immediately.
- 2) For Naïve Cells with Short Stim, Take 10 million and stimulate in T cell media at  $37^\circ$  and 7%  $\text{CO}_2$  with plate-bound anti-CD3 + anti-CD28 for 4 hours. Then freeze cells.
- 3) For Th0 Cells, take 10 million and culture with CD3/CD28 beads with no cytokines added (neutral conditions) in T cell media at  $37^\circ$  and 7%  $\text{CO}_2$  for 3 days. Remove beads, collect cells, and freeze.
- 4) For Th1 Cells, take 10 million and culture with CD3/CD28 beads + 1 ng/ml IL-12 (Th1 conditions) in T cell media at  $37^\circ$  degrees and 7%  $\text{CO}_2$  for 3 days. Remove beads, collect cells, and freeze.